


**RABBIT ANTI- CYP450 1A1
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB10088 **QUANTITY:** 50 µL
LOT NUMBER: **HOST/ISOTYPE:** Rabbit
IMMUNOGEN: Recombinant CYP450 2A1
APPLICATIONS: Western Blot Analysis: Western Blot Analysis of AB10080
at a dilution of 1:5,000.

Lane 1 Recombinant CYP450 1A1
Lane 2 Napthalone induced Rat Microsomes
Arrow indicates expected band at 58 kDa



Optimal working dilutions must be determined by the end user.

SPECIES REACTIVITY: Human
FORMAT: Raw Rabbit Serum
PRESENTATION: Raw Rabbit Serum with 0.1% Sodium Azide
STORAGE/HANDLING: Maintain at 2-8°C for up to 12 months from date of receipt.
REFERENCES: Thapliyal, R and Maru, G B (2001). Inhibition of cytochrome P450 isozymes by curcumins in vitro and in vivo. *Food Chem Toxicol* **39**: 541-7.
Smith, Cornelia M, *et al.* Modulation of UDP-glucuronosyltransferase 1A1 in primary human hepatocytes by prototypical inducers. *J Biochem Mol Toxicol* **19**: 96-108.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.