

**ANTI-HCV CORE
MONOCLONAL ANTIBODY**

CATALOG NUMBER:	MAB8682	LOT NUMBER:	
QUANTITY:	100 µg	CONCENTRATION:	0.75 mg/mL
CLONE NAME:	1F6	HOST/ISOTYPE:	Mouse/IgG _{2a}
SPECIFICITY:	HCV core antigen. Recognizes 80-120 amino acid residues of HCV core antigen		
IMMUNOGEN:	Produced against a recombinant protein		
SPECIES REACTIVITY:	Reacts with human. Reactivity with other species has not been confirmed.		
APPLICATIONS:	ELISA: 1:1000 Does not work for Western Blotting Optimal working dilutions must be determined by end user.		
FORMAT:	Affinity purified		
PRESENTATION:	Liquid in phosphate buffered saline, pH 7.2 with 0.01% sodium azide.		
STORAGE/HANDLING:	Store at -20°C in undiluted aliquots. May be stored at 2° to 8°C for short-term use. Avoid repeated freeze-thaw cycles.		

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.