

Get the Access Services Advantage

Our dedicated team has the in-depth product knowledge and the application expertise to deliver quality results every time, right on schedule. With 15 years of

experience and research to draw on, there aren't many microbiological testing or validation challenges we haven't seen – and solved.

In addition to the time and cost savings, you'll appreciate the convenience and confidence that come with Millipore's Access Services.



Milliflex™ Test System

▶ Bioburden Testing

▶ Microbial Limit Testing

▶ Environmental Monitoring

▶ Sterility Testing

MicropreSure™ Filter Units

M Air T™ Air Monitor

Steritest™ Disposable Devices

Steritest® Sterility Testing System

To discuss a specific application, please contact your Millipore representative, or call **1-800-MILLIPORE**.

To Place an Order or Receive Technical Assistance

For additional information, call your nearest Millipore office.

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

Internet: www.millipore.com

Tech Service: www.millipore.com/techservice



Now you can buy Millipore products online @ www.millipore.com/purecommerce



For more information about Millipore's Access Quality Assurance Services, or about any of our pharmaceutical Quality Assurance and Quality Control solutions, contact the Millipore subsidiary nearest you, or call 781-533-6000.

In the U.S. and Canada, call toll-free **1-800-MILLIPORE**.

You'll also find a wealth of information on our web site at www.millipore.com

MILLIPORE

Millipore and Steritest are registered trademarks of Millipore Corporation. Milliflex, MicroStar, M Air T and Sterisure are trademarks of Millipore Corporation. Access is a service mark of Millipore Corporation. PB1150EN00 Printed in U.S.A. 02/02 01-179 Copyright 2002 Millipore Corporation Bedford, MA 01730 U.S.A. All rights reserved.

Millipore Worldwide

- AUSTRALIA**
Tel. 1 800 222 111
or (02) 9868 8999
Fax (02) 9878 0788
- AUSTRIA**
Tel. (01) 877 89 26
Fax (01) 877 16 54
- BALTIC COUNTRIES**
Tel. +358 9 804 5110
Fax +358 9 256 5660
- BELGIUM AND LUXEMBOURG**
Tel. +32 2 226 88 40
Fax +32 2 726 98 84
- BRAZIL**
Tel. (011) 5548-7011
Fax (011) 5548-7923
- CANADA**
Tel. 1-800-645-5476
Fax 1-800-645-5439
- CHINA, PEOPLE'S REPUBLIC OF**
Beijing:
Tel. (8610) 8519 1250
(8610) 6518 1058
Fax (8610) 8519 1255
Guangzhou:
Tel. (8620) 8755 4049
Fax (8620) 8752 0172
Hong Kong:
Tel. (852) 2803 9111
Fax (852) 2513 0313
Shanghai:
Tel. (8621) 5306 9100
Fax (8621) 5306 0838
- CZECH REPUBLIC**
Tel. 02-2051 3841
Fax 02-2051 4298
- DENMARK**
Tel. 70 10 00 23
Fax 70 10 13 14
- EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF**
Life Sciences Division:
Tel. +33 3 88 38 9536
Fax +33 3 88 38 9539
BioPharmaceutical Division:
Tel. +43 1 877-8926
Fax +43 1 877-1654
Laboratory Water Division:
Tel. +33 1 30 12 7000
Fax +33 1 30 12 7180
- FINLAND**
Tel. (09) 804 5110
Fax (09) 256 5660
- FRANCE**
Tel. 01 30 12 7000
Fax 01 30 12 7180
- GERMANY**
Tel. (06196) 494-0
Fax (06196) 43901
- HUNGARY**
Tel. 01 381 0433
01 381 0434
01 209 3232
Fax 01 209 0295
- INDIA**
Tel. (91) 80-839 46 57
Fax (91) 80-839 63 45
- IRELAND**
Life Sciences and Laboratory Water Divisions:
Tel. +44 1923 816375
Fax +44 1923 818297
BioPharmaceutical Division:
Tel. (021) 883 666
Fax (021) 883 048
- ITALY**
Vimodrone (Milano):
Tel. (02) 25 07 81
Fax (02) 26 50 324
Roma:
Tel. (06) 52 03 600
Fax (06) 52 95 735
- JAPAN**
Tel. (03) 5442-9711
Fax (03) 5442-9736
Life Sciences 9737
BioPharm. 9734
Lab Water
- KOREA**
Tel. (822) 551-0990
Fax (822) 551-0228
- MALAYSIA**
Tel. 603-7957-1322
Fax 603-7957-1711
- MEXICO**
Tel. (525) 576-9688
Fax (525) 576-8706
- THE NETHERLANDS**
Tel. 076-5022000
Fax 076-5022436
- NORWAY**
Tel. 22 67 82 53
Fax 22 66 04 60
- POLAND**
Tel. 22-669 12 25
22-663 70 31
Fax 22-663 70 33
- PUERTO RICO**
Tel. (787) 273-8495
Fax (787) 747-6553
- SINGAPORE**
Tel. (65) 842 1822
Fax (65) 842 4988
- SPAIN AND PORTUGAL**
Tel. +34 917 283 900
Fax +34 917 292 909
- SWEDEN**
Tel. 08-628 6960
Fax 08-628 6457
- SWITZERLAND**
Tel. (01) 908-30-60
Fax (01) 908-30-80
- TAIWAN**
Tel. 886-2-2792-9333
Fax 886-2-2792-6555
- U.K.**
Tel. 01923 816375
Fax 01923 818297
- U.S.A.**
Tel. (781) 533-6000
Fax (781) 533-3110
- IN ALL OTHER COUNTRIES**
Millipore Intertech (U.S.A.)
Tel. +1 (781) 533-8622
Fax +1 (781) 533-8630



pure assurance

Access™ QA/QC Services

Cost savings, convenience, compliance and confidence

Quality services complement our quality products



A complete quality assurance solution means more than just reliable products. It means getting your testing methods and equipment validated, and producing results as quickly and accurately as possible.

Millipore's Access Services was established in 1987 to give resource-constrained QA/QC groups a welcome alternative to in-house methods development and validation. Whether you need a new method or protocol, on-site IQ, OQ or PQ, or a complete, company-customized validation package, we're ready to lend a hand.

Save time and money

We can help develop streamlined methods that save you steps, time and money. And we'll significantly shorten your time to validation with professional, pre-packaged validation protocols that include everything you need to achieve full compliance.

One less item for your to-do list

By taking care of the details needed to get your methods and equipment validated and fully operational, we can free up your staff for other important projects.

Compliance – and confidence

Millipore Validation Packages and expertise can reduce your risk, and make the job much faster and easier. With years of experience helping pharmaceutical and biotech companies successfully meet regulatory requirements, we're familiar with current guidelines from the FDA, EMEA, and other regulatory authorities.

We not only share your sense of urgency, we share your insistence on getting the job done right. With Millipore, you'll know – for sure – that you're getting top performance from your equipment, in full regulatory compliance.

Laboratory Services and Methods Development

Millipore test methods are used in QA/QC labs of many leading pharmaceutical and biotech manufacturers around the globe. If you need to develop a new method or optimize an existing protocol, we can help you do it quickly and easily – in strictest confidence, of course.

Our scientists are available for consultation and project-based laboratory work at our state-of-the-art facilities. Typical assignments include:

- Feasibility and optimization studies
- Development of microbiological test methods to assay a variety of biopharmaceuticals
- Neutralization of antimicrobial additives or procedures
- Product compatibility testing

Standard Validation Packages

Millipore's standard Validation Packages are an economical way to meet the needs of most laboratories, and are ideal for busy QA/QC departments. They not only save you money; they can lighten your workload, reduce paperwork, and get you up and running faster.

For many of our QC testing solutions – such as Milliflex, Steritest, M Air T, Sterisure and MicroStar systems – we provide a complete, ready-to-run validation package that consists of:

- A validation master plan
- A detailed specification
- Standard IQ and OQ protocols
- PQ guidelines that can be easily customized for your specific application
- A Validation Final Report format, with all necessary worksheets



Customized Validation Packages

Our validation experts will work closely with your validation team to create customized validation packages to meet highly individual needs. These are based on our Standard Validation Package, reformatted to conform to your company's specifications, or to incorporate required modifications to test methods.

On-Site Validation

To save time, and minimize disruption of your normal workflow, we can also help you execute your validation study. All work is performed at your location by an experienced Validation Specialist. With validation complete, your test equipment will be ready for use.

▶ *Lighten your workload*

▶ *Save time and money*

▶ *Get up and running faster*

