


RABBIT ANTI-GRASP-1 AFFINITY PURIFIED POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB9636	QUANTITY:	100 µL
LOT NUMBER:			
ALTERNATE NAMES:	GRIPAP1	EPITOPE:	C-Terminal

SPECIFICITY: GRASP-1 (General receptor for phosphoinositides 1-associated scaffold protein, GRP1-associated scaffold protein). The antibody recognizes a protein of ~95 kDa corresponding to GRASP-1 in lysates from rat brain cerebellum.

IMMUNOGEN: Synthetic peptide from the C-terminal of rat GRASP-1.

APPLICATIONS: Western blot: 1:1,000



Optimal working dilutions must be determined by the end user.

SPECIES REACTIVITY: Rat. Other species have not been tested. The immunogen has 100% homology with human and mouse.

FORMAT: Affinity purified immunoglobulin

PRESENTATION: Liquid in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg/mL BSA and 50% glycerol.

STORAGE/HANDLING: Maintain at -20°C in undiluted for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles. Do not store in a self defrosting freezer.

REFERENCES: Scannevin, R.H., Haganir, R.L., *Nat. Rev. Neurosci.*(2000) 1:133-141.
Ye, B., et al., *Neuron* (2000) 26:603-617.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com