

**RABBIT ANTI-SUBSTANCE P
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB1977

LOT NUMBER:

QUANTITY: 50 µL

SPECIFICITY: Substance P.

IMMUNOGEN: Substance P, BSA conjugate.

APPLICATIONS: Immunohistochemistry* on human intestine (ABC): 1:5,000-1:10,000.
* fixative – Zamboni's
Optimal working dilutions must be determined by the end user.

INHIBITION CONTROL: 10 nmol of Substance P per mL of diluted antiserum.

FORMAT: Neat rabbit serum

PRESENTATION: Liquid, no preservatives.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Have you published a paper using a Chemicon primary antibody? You may be eligible for credit toward future purchases of Chemicon products!
Contact our Technical Service department for more information.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.