



ANTI-HUMAN IgA,IgG,IgM, (H&L chains)
IgG FRACTION, FITC CONJUGATED, POLYCLONAL ANTIBODY
CERTIFICATE OF ANALYSIS

CATALOG NUMBER:	AB7114F	QUANTITY:	2 mL
LOT NUMBER:		CONCENTRATION:	
HOST/ISOTYPE:	Sheep, IgG		
SPECIFICITY:	Prior to conjugation, the antibody was specific for immunoglobulins when tested against human plasma and a range of monoclonal immunoglobulins in pathological plasma and urine samples by Coomassie Blue R stained double immunodiffusion and immunoelectrophoresis and found to be specific for human immunoglobulins.		
IMMUNOGEN:	Purified human immunoglobulins from pooled human serum.		
APPLICATIONS:	IF: 1/10-1/400 <i>Optimal dilutions to be determined by enduser and exact assay.</i>		
SPECIES REACTIVITY:	Human, other species not examined.		
PRESENTATION:	Liquid in 20 mM Tris, 70 mM NaCL pH 8.0, with 0.1% sodium azide as a preservative..		
STORAGE/HANDLING:	Store at 2-8°C protected from light and temperature fluxuations for up to expiration date given on product label.		

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals
©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com