

RABBIT ANTI-OVALBUMIN (Hen Egg White) POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1225
LOT NUMBER:	
QUANTITY:	25 mg
CONCENTRATION:	10 mg/mL
SPECIFICITY:	Ovalbumin (Hen Egg White). Monospecific by IEP.
IMMUNOGEN:	Ovalbumin (Hen Egg White)
APPLICATIONS:	Western blot: 1:500 – 1:2000 Immunohistochemistry: 1:300-1:2,000 ELISA: 1:3,000-1:12,000 when assayed against 1 µg of Ovalbumin in a standard sandwich Optimal working dilutions must be determined by end user.
FORMAT:	Immunoglobulin purified by delipidation, salt fractionation and ion-exchange chromatography.
PRESENTATION:	Liquid in 0.02 M Potassium Phosphate buffer, 0.15 M Sodium Chloride, pH 7.2, with 0.01% sodium azide.
STORAGE/HANDLING:	Maintain frozen at -20°C in undiluted aliquots for up to six months after date of receipt. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003-2006: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com

2-10-2006/AB1225/CB
/10/24/07/ft