



CHEMICON IHC SELECT RESEARCH®
PREDILUTED RABBIT ANTI-TYROSINE HYDROXYLASE
AFFINITY PURIFIED
POLYCLONAL ANTIBODY

CATALOG NUMBER:	IHCR1003-6	QUANTITY:	6 mL
LOT NUMBER:	TR1443115		
DESCRIPTION:	Tyrosine hydroxylase plays an important role in the physiology of adrenergic neurons. It is the first rate-limiting enzyme involved in the biosynthesis of the catecholamines Dopamine and Norepinephrine from tyrosine. TH is, therefore, a useful marker for dopaminergic and noradrenergic neurons. The enzymatic activity of TH requires ferrous ions as cofactors and is believed to be regulated by phosphorylation. At least four isoforms of human TH have been identified which results from alternative splicing.		
SPECIFICITY:	Tyrosine Hydroxylase		
IMMUNOGEN:	Denatured tyrosine hydroxylase from rat pheochromocytoma (denatured by sodium dodecyl sulfate).		
APPLICATIONS:	Antibody is prediluted and ready to use for Immunohistochemistry of formalin-fixed, paraffin-embedded tissues. <u>Pretreatment:</u> Heat Induced Epitope Retrieval (HIER). Recommend Citrate Buffer, pH 6.0 (Cat. No. 21545) or EDTA, pH 8.0 (Cat. No. 21546). No signal was detected without Epitope retrieval. <u>Incubation:</u> 30 minutes with IHC Select® Detection Kits. Tyrosine Hydroxylase has been prediluted for use as the primary antibody with Chemicon's IHC Select® Detection Kits and Protocols (Catalog Nos. DAB050, DET-HP1000, APR050, and DET-APR1000), but other supplier's IHC detection systems may be used. For optimized protocol details, visit www.chemicon.com and select the protocols link under Cat. No.IHCR1003-6.		
SPECIES REACTIVITIES:	Reacts with Human. Reactivity with other species has not been determined.		
PRESENTATION:	Liquid diluted in PBS, pH 7.2 with stabilizers, 0.2% Tween 20, and 0.1% ProClin 300 as preservative.		
STORAGE/HANDLING:	Maintain at 2-8°C. Refer to vial for expiration dating.		

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.



**RELATED
REFERENCES:**

- Hollis, J. et al. (2005). Lipopolysaccharide has selective actions on sub-populations of catecholaminergic neurons involved in activation of the hypothalamic–pituitary–adrenal axis and inhibition of prolactin secretion. *J. Endocrinology* **184** : 393-406.
- Witkovsky, P. et al. (2004). Activity-dependent phosphorylation of tyrosine hydroxylase in dopaminergic neurons of the rat retina. *J. Neurosci.* **24(17)** : 4242-4249.
- Alberi, L et al. (2004). Engrailed genes are cell-autonomously required to prevent apoptosis in mesencephalic dopaminergic neurons. *Development* **131(13)** : 3229-3236.
- Baker, SA et al. (2004). Dopaminergic nigrostriatal projections regulate neural precursor proliferation in the adult mouse subventricular zone. *Eur J Neurosci* **20(2)** : 575.
- Murray, TK et al. (2003). LY503430, a novel AMPA receptor potentiator with functional, neuroprotective and neurotrophic effects in rodent models of Parkinson's disease.. *A Pharma Exp Therap JPET* **49445**
- Contini, M. and Raviola, E. (2003). GABAergic synapses made by a retinal dopaminergic neuron . *Proc. Nat. Acad. Sci. (USA)* **100(3)** : 1358-1363.
- Wu, P. et al. (2002). Region-specific generation of cholinergic neurons from fetal human neural stem cells grafted in adult rat. *Nature Neurosci* **5(12)** : 1271-1278.
- Tillerson, JL. et al. (2002). Forced Nonuse in Unilateral Parkinsonian Rats Exacerbates Injury. *J Neurosci* **22(15)** : 6790-6799.
- Prohaska, J. and Brokate, B. (2001). Dietary copper deficiency alters protein levels of rat dopamine beta-monooxygenase and tyrosine monooxygenase. *Exp. Biol. Med.* **226(3)** : 199-207.
- Feany, M. and Bender, W. (2000). A Drosophila model of Parkinson's disease. *Nature* **404** : 394-398.
- Mirabella, N., et al. (2000). The distribution and co-localization of nicotine adenine dinucleotide phosphate reduced-diaphorase (NAH-d)-/nitric oxide synthase (NOS)-containing neurons in the innervation of the duck ureter. *Anat.Embryol. (Berl)* **202** : 291-301.
- Lode, H., et al. (2000). Tyrosine hydroxylase-based DNA-vaccination is effective against murine neuroblastoma. *Medical and Pediatric Oncology* **35** : 641-646.
- Nichol, K., et al. (1999). Characterization of hypothalamic neurons expressing a neuropeptide receptor, GALR2, using combined *in situ* hybridization-immunohistochemistry. *Methods* **18** : 481-486.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

- Fernandez-Galaz, M., et al. (1999). Diurnal oscillation in glial fibrillary acidic protein in a perisuprachiasmatic area and its relationship to the luteinizing hormone surge in the female rat. *Neuroendocrinology* **70** : 368-376.
- Jaber, M., et al. (1999). Differential regulation of tyrosine hydroxylase in the basal ganglia of mice lacking the dopamine transporter. *Eur J. Neurosci.* **11(10)** : 3499-3511.
- Lew, J., et al. (1999). Increased site-specific phosphorylation of tyrosine hydroxylase accompanies stimulation of enzymatic activity induced by cessation of dopamine neuronal activity. *Mol. Pharm.* **55(2)** : 202-209.
- Shinder, V et al. (1999). Structural basis of sympathetic-sensory coupling in rat and human dorsal root ganglia following peripheral nerve injury. *J Neurocytology* **28(9)** : 743-761.
- Horger, B., et al. (1998). Neurturin exerts potent actions on survival and function of midbrain dopaminergic neurons. *Journal of Neuroscience* **18(13)** : 4929-4937.
- Srivatsan, M. ,and Peretz, B. (1997). Acetylcholinesterase promotes regeneration of neurites in cultured adult neurons of *Aplysia*. *Neuroscience* **77(3)** : 921-931.
- Snyder-Keller, A. and Costantini, L. (1996). Glutamate receptor subtypes localize to patches in the developing striatum. *Developmental Brain Research* **94** : 246-250.
- Giordano, A., et al. (1996). Tyrosine hydroxylase, neuropeptide Y, substance P, calcitonin gene-related peptide and vasoactive intestinal peptide in nerves of rat periovarian adipose tissue: an immunohistochemical and ultrastructural investigation. *Journal of Neurocytology* **25** : 125-136.
- Davis, B., et al. (1996). Effects of NGF overexpression on anatomical and physiological properties of sympathetic postganglionic neurons. *Brain Research* **724** : 47-54.
- Carlson, S., et al. (1995). NGF modulates sympathetic innervation of lymphoid tissues. *Journal of Neuroscience* **15(9)** : 5892-5899.
- Rayport, S. and Sulzer, D. (1995). Visualization of antipsychotic drug binding to living mesolimbic neurons reveals D2 receptor, acidotropic, and lipophilic components. *Journal of Neurochemistry* **186** : 691-703.
- Harper, M. and Wallace, M. (1995). Changes in density of brainstem afferents in ferret primary auditory cortex (AI) during postnatal development. *Journal of Anatomy* **186** : 373-382.
- Albers, K., et al. (1994). Overexpression of nerve growth factor in epidermis of transgenic mice causes hypertrophy of the peripheral nervous system. *Journal of Neuroscience* **14(3)**

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209
www.millipore.com



: 1422-1432.

Barnett, E., et al. (1993). Two neurotrophic viruses, herpes simplex virus type 1 and mouse hepatitis virus, spread along different neural pathways from the main olfactory bulb. *Neuroscience* **57(4)** : 1007-1025.

Ren, T. and Sagar (1992). Induction of c-fos immunostaining in the rat brain after the systemic administration of nicotine. *Brain Research Bulletin* **29** : 589-597.

Wang, Z.-Z., et al. (1991). Co-existence of tyrosine hydroxylase and dopamine beta-hydroxylase immunoreactivity in glomus cells of the cat carotid body. *Journal of the Autonomic Nervous System* **32** : 259-264.

Haycock, J. W. (1989). Quantitation of tyrosine hydroxylase protein levels: spot immunolabeling with an affinity-purified antibody. *Anal. Biochem.* **181(2)** : 259-266.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.