

**RABBIT ANTI-UTF1
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB3383

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SPECIFICITY: UTF1

IMMUNOGEN: Synthetic peptide from human UTF1.

APPLICATIONS: Western blot: 0.5-1.0 µg/mL on mouse ES cell extract. Reacts with a band of ~36 kDa. Also reacts with a band of ~23 kDa. Optimal working dilutions must be determined by the end user.

SPECIES REACTIVITIES: Mouse. Other species have not been tested.

FORMAT: Affinity purified immunoglobulin.

PRESENTATION: Liquid in 0.02M phosphate buffer containing 0.25M NaCl, pH 7.6 with 0.1% sodium azide.

STORAGE/HANDLING: Maintain at 2-8°C in undiluted for up to 6 months after date of receipt.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com