

**RABBIT ANTI-MOUSE VASCULAR ENDOTHELIAL GROWTH FACTOR (VEGF)
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB1876

LOT NUMBER: xxx

QUANTITY: 100 µg

CONCENTRATION: 0.90 mg/mL

SPECIFICITY: Reacts with mouse VEGF. No cross-reactivity with WHO standards to an extensive panel of cytokines and growth factors by EIA.

IMMUNOGEN: Recombinant mouse VEGF.

APPLICATIONS: Western blot: 1:100-1:1,000
ELISA: 1:100-1:1,000
Optimal working dilutions must be determined by the end user.

FORMAT: Immunoglobulin purified by HPLC.

PRESENTATION: Liquid (sterile filtered) with no preservatives.

STORAGE/HANDLING: Maintain at 2-8°C in undiluted aliquots for up to 12 months. Dilute with Tris Buffered Saline (TBS).

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial on a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com