

**RABBIT ANTI- DOPAMINE RECEPTOR D5  
AFFINITY PURIFIED  
POLYCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	AB9509
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	50 µg
<b>CONCENTRATION:</b>	1.0 mg/mL
<b>SPECIFICITY:</b>	Recognizes Dopamine Receptor D5.
<b>IMMUNOGEN:</b>	Synthetic peptide from the 3rd extracellular domain of human Dopamine Receptor D5.
<b>APPLICATIONS:</b>	Immunohistochemistry on formalin fixed paraffin embedded tissue sections. Optimal working dilutions must be determined by the end user.
<b>SPECIES REACTIVITIES:</b>	Human and rat. Other species have not been tested. The immunogen has 62% homology with mouse and 62% with rat.
<b>FORMAT:</b>	Affinity purified immunoglobulin.
<b>PRESENTATION:</b>	Liquid in PBS, pH 7.7 with 0.1% sodium azide.
<b>STORAGE/HANDLING:</b>	Maintain at -70°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

---

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000  
www.millipore.com