

**RABBIT ANTI-BRAIN LIPID-BINDING PROTEIN (BLBP)
POLYCLONAL ANTIBODY**

- CATALOG NUMBER:** AB9558
- LOT NUMBER:**
- QUANTITY:** 100 μ L
- SPECIFICITY:** Recognizes Brain Lipid-Binding Protein (BLBP). The antibody is specific to radial glia (embryonically) and astrocytes postnatally.
- IMMUNOGEN:** Recombinant whole BLBP.
- APPLICATIONS:** Immunoblotting: 1:1,000-1:5,000 using ECL. Reacts with the ~15 kD BLBP. Immunohistochemistry: 1:500-1:1,500. The antibody works on cryostat/vibratome sectioned tissue. It has been reported that the antibody works on paraffin-embedded materials. Optimal working dilutions must be determined by the end user.
- SPECIES REACTIVITIES:** Rat and mouse. Other species have not been tested.
- FORMAT:** Rabbit serum.
- PRESENTATION:** Liquid. Contains no preservative.
- STORAGE/HANDLING:** Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
- REFERENCE:** Feng, L., et al., *Neuron* (1994) **12**:895-908.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.