

HISTAMINE RECEPTOR 1 CONTROL PEPTIDE

- CATALOG NUMBER:** AG204
- LOT NUMBER:**
- QUANTITY:** 100 µg
- CONCENTRATION:** 1 mg/mL
- DESCRIPTION:** A 19 amino acid peptide sequence within the 3rd cytoplasmic loop #3 of human HR1 (1).
- APPLICATIONS:** Antibody blocking (AB5654P)
Optimal working dilution must be determined by the end user.
- PRESENTATION:** Liquid in PBS containing 0.05% sodium azide.
- STORAGE/HANDLING:** Store at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
- RELATED REFERENCES:** 1) Fujimoto K et al (1993) *BBRC*, **190**:294; Inoue I et al (1996) *PNAS* **93**:13316; DeBacker MD et al (1993) *BBRC* **197**:1601; Fukui H et al (1994) *BBRC* **201**:894; Moguilevsky N et al (1994) *Eur. J. Biochem.* **224**:489.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

© 2003-2004 - CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.