

**MOUSE ANTI-SV40 T ANTIGEN
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB986

LOT NUMBER:

QUANTITY: 100 µL

SPECIFICITY: SV40 T antigen. Precipitates large T antigen and NVT antigen (when complexed) but not small t antigen (nor SV-80-specific smaller t antigen) from extracts of SV-80 cells and of SV-40 infected TC-7 cells. Recognizes denaturation-sensitive antigenic determinants unique to large T antigen. Does not bind gel-purified denatured large T antigen or mouse NVT antigen in RIA.

IMMUNOGEN: SV40 transformed mouse cells.

ISOTYPE: IgG₁

APPLICATIONS: Dilute with buffer pH 7.5 - 8.0 to desired working volumes. For extensive dilution, protein-containing or other stabilizing medium should be used. Final working dilutions must be determined by end user.

FORMAT: Ascites.

STORAGE/HANDLING: Maintain at -20°C for up to 12 months. Avoid repeated freeze/thaw cycles.

REFERENCES: *J. Virology* **34**:752-763 (1980).

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com