

**METABOTROPIC GLUTAMATE RECEPTOR 2/3 (mGluR2/3)
CONTROL PEPTIDE**

CATALOG NUMBER:	AG212
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1.0 mg/mL
DESCRIPTION:	Peptide from the carboxy terminus of rat mGluR2
APPLICATIONS:	For use in blocking the reactivity of AB1553. Optimal working dilution must be determined by the end user.
PRESENTATION:	Liquid in water. Contains no preservative.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.