

GOAT ANTI-APOLIPOPROTEIN B
POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB742
LOT NUMBER:	
QUANTITY:	1 mL
SPECIFICITY:	Apolipoprotein B. Monospecific by IEP with neat antibody vs. 2X pooled human serum and neat pooled human plasma.
IMMUNOGEN:	Purified human Apolipoprotein B.
FORMAT:	Defibrinated, delipidized and absorbed by solid phase chromatography as required.
PRESENTATION:	Goat antisera in 0.05M Tris-HCl pH 7.5, 0.5M NaCl with 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCES:	Liao, W., et al., <i>Biochemistry</i> (1999) 38 :7532-7544. Liao, W., et al., <i>Journal of Biological Chemistry</i> (1998) 273 :27225-27230. Rashid, K.A., et al., <i>Journal of Biological Chemistry</i> (2002) 277 :22010-22017.

For research use only; not for use as a diagnostic.