

**MOUSE ANTI-SYNAPTOBREVIN [VAMP]
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB335

LOT NUMBER:

VOLUME: 100 μ L

SPECIFICITY: Synaptobrevin [VAMP], 16kDa protein. Not reactive with bacterially expressed fusion protein.

IMMUNOGEN: Crude synaptic immunoprecipitate (Human)

ISOTYPE: IgG₁

APPLICATIONS: Immunohistochemistry at 1:100-1:200
Formalin fixed (fixed <4 weeks), paraffin embedded
Frozen
4% Paraformaldehyde
Western blotting at 1:100-1:1,000.
ELISA at 1:500
Optimal working dilutions must be determined by the end user.

SPECIES REACTIVITIES: Human, cat and pig.

FORMAT: Purified immunoglobulin.

PRESENTATION: Liquid in 0.02M Phosphate buffer, 0.25M NaCl with 0.1% sodium azide.

STORAGE/HANDLING: Maintain at 2-8°C in undiluted aliquots for up to 6 months after date of receipt.

REFERENCES: *Brain Research* **609**:9-20 (1993).
Acta Neuropathol **87**:135-142 (1994).

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial on a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.