

CHEMISCREEN™ MEMBRANE PREPARATION RECOMBINANT HUMAN GLUCAGON RECEPTOR

CATALOG NUMBER:	HTS112M	QUANTITY:	200 units
LOT NUMBER:		VOLUME/CONCENTRATION :	1 mL, 1 mg/mL

BACKGROUND: Glucagon is a 29-amino acid peptide that stimulates glycogenolysis and gluconeogenesis in the liver to increase blood glucose. The receptor for glucagon is a class 2 (or class B) GPCR that signals through G_s to stimulate cAMP production (Mayo *et al.*, 2003). Mice lacking the glucagon receptor have mild hypoglycemia after fasting, and exhibit hyperplasia of pancreatic α-cells (Gelling *et al.*, 2003). Because of its role in promoting hyperglycemia, the glucagon receptor presents a potential target for treatment of diabetes. Millipore's glucagon receptor membrane preparations are crude membrane preparations made from our proprietary stable recombinant cell lines to ensure high-level of GPCR surface expression; thus, they are ideal HTS tools for screening of antagonists of glucagon receptor interactions and its ligands. The membrane preparations exhibit a K_d of 3.09 nM for [¹²⁵I]-glucagon. With 1 nM [¹²⁵I]-glucagon, 5 μg/well Glucagon Receptor Membrane Prep yields greater than 10-fold signal-to-background ratio.

APPLICATIONS: Radioligand binding assay




Figure 1. Saturation binding for Glucagon. 5 μg/well Glucagon Receptor Membrane Preparation was incubated with increasing amount of ¹²⁵I labeled glucagon in the absence (total binding, TB) or presence (nonspecific binding, NSB) of 500-fold excess unlabeled glucagon. Specific binding (SB) was determined by subtracting NSB from TB.




Figure 2. Competition binding for Glucagon. Glucagon Receptor Membrane Preparation (5 and 10 μg/well in a 96-well plate) were incubated with 1 nM ¹²⁵I labeled glucagon and increasing concentrations of unlabeled glucagon, and subjected to filtration binding.

Table 1. Signal:background and specific binding values obtained in a competition binding assay with varying amounts of Glucagon receptor membrane prep.

	10 µg/well	5 µg/well
Signal:background	10.1	35
Specific binding (cpm)	33315	36252

SPECIFICATIONS: 1 unit = 5 µg

 B_{max} for [¹²⁵I] glucagon binding: 16.7 pmol/mg protein; K_d for [¹²⁵I] glucagon binding: ~ 3.09 nM

TRANSFECTION: Full-length human GCGR cDNA encoding Glucagon Receptor (Accession Number: NM_000160)

HOST CELLS: Chem-1, an adherent mammalian cell line with no endogenous glucagon receptor expression.

RECOMMENDED ASSAY CONDITIONS: Membranes are mixed with radioactive ligand and unlabeled competitor (see Figures 1 and 2 for concentrations tested) in binding buffer in a nonbinding 96-well plate, and incubated for 1-2 h. Prior to filtration, an FC 96-well harvest plate (Millipore cat. # MAHF C1H) is coated with 0.33% polyethyleneimine for 30 min, then washed with 50 mM HEPES, pH 7.4, 0.5% BSA. Binding reaction is transferred to the filter plate, and washed 3 times (1 mL per well per wash) with Wash Buffer. The plate is dried and counted.

Binding buffer: 50 mM Hepes, pH 7.4, 5 mM MgCl₂, 1 mM CaCl₂, 0.2% BSA, filtered and stored at 4°CRadioligand: [¹²⁵I]-glucagon (Perkin Elmer#:NEX-207)

Wash Buffer: 50 mM Hepes, pH 7.4, 500mM NaCl , 0.1% BSA, filtered and stored at 4°C .

One package contains enough membranes for at least 200 assays (units), where a unit is the amount of membrane that will yield greater than 10-fold signal:background with ¹²⁵I labeled glucagon at 1 nM**PRESENTATION:**

Liquid in packaging buffer: 50 mM Tris pH 7.4, 10% glycerol and 1% BSA with no preservatives.

Packaging method: Membranes protein were adjusted to the indicated concentration in packaging buffer, rapidly frozen, and stored at -80°C.

STORAGE/HANDLING:

Store at -70°C. Product is stable for at least 6 months from the date of receipt when stored as directed. Do not freeze and thaw.

REFERENCES:Gelling RW *et al.* (2003) Lower blood glucose, hyperglucagonemia, and pancreatic α cell hyperplasia in glucagon receptor knockout mice. *Proc. Natl. Aca. Sci. USA* 100:1438-1443.Mayo KE *et al.* (2003) International Union of Pharmacology. XXXV. The glucagon receptor family. *Pharmacol. Rev.* 55: 167-194.*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2006-2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com