

**MOUSE ANTI- HUMAN METAPNEUMOVIRUS (MPV)
MONOCLONAL ANTIBODY**

CATALOG NUMBER:	MAB8510		
LOT NUMBER:			
QUANTITY:	100 µg	CONCENTRATION:	1 mg/mL
CLONE NAME:	1B7	HOST/ISOTYPE:	IgG ₂
SPECIFICITY:	Reacts against both A and B antigenic groups and subgroups 1 and 2		
IMMUNOGEN:	MPV CAN98-75		
APPLICATIONS:	ELISA IF Optimal dilutions must be determined by the end user		
SPECIES REACTIVITY:	Human Metapneumoviruses		
FORMAT:	Purified immunoglobulin.		
PRESENTATION:	Liquid in 0.02 M Phosphate buffer, 0.25 M Sodium Chloride, 0.1% sodium azide, pH = 7.6.		
STORAGE/HANDLING:	Store at 2° to 8°C for up to 12 months after date of receipt		
REFERENCES:	Bastien, N. et al. 2003. Sequence analysis of the N, P, M and F genes of Canadian human metapneumovirus strains. <i>Virus Res.</i> May <u>93</u> (1):51-62 Boivin, G. et al. 2002. Virological features and clinical manifestations associated with human metapneumovirus: a new paramyxovirus responsible for acute respiratory-tract infections in all age groups. <i>J. Infect. Dis.</i> <u>186</u> : 1330-1334.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com