

HISTAMINE RECEPTOR 4 (HR4) CONTROL PEPTIDE

CATALOG NUMBER:	AG217
LOT NUMBER:	22030700
QUANTITY:	100 µg
CONCENTRATION:	1.0 mg/mL
DESCRIPTION:	Synthetic 17 amino acid peptide sequence within the cytoplasmic, C-terminus of human HR4.
APPLICATIONS:	For use in blocking the reactivity of AB5663P. Optimal working dilution must be determined by the end user.
PRESENTATION:	Liquid in PBS containing 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles
RELATED REFERENCES:	Oda T et al (2000) <i>J. Biol. Chem.</i> 275 :36781-36786; Nakamura T et al (2000) <i>BBRC</i> 279 :615-620. <i>For research use only; not for use as a diagnostic.</i>
Important Note:	<i>During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.</i>