

**RABBIT ANTI-BOVINE SURFACTANT PROTEIN A (SP-A)
POLYCLONAL ANTIBODY**

CATALOG NO.: AB3424

LOT NUMBER:

QUANTITY: 100 µL

SPECIFICITY: Reacts strongly with bovine and sheep SP-A.

IMMUNOGEN: Bovine SP-A

APPLICATIONS: Immunoblot: 1:1,000-1:5,000
ELISA: 1:1,000-1:30,000
Optimal dilutions must be determined by the end user.

FORMAT: Rabbit serum.

PRESENTATION: Liquid. Contains no preservative.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com