

## FIXATION SOLUTION, 5X

---

<b>CATALOG NUMBER:</b>	5113	<b>QUANTITY:</b>	220 mL
<b>LOT NUMBER:</b>		<b>CONCENTRATION:</b>	5X
<b>APPLICATIONS:</b>	For use with Catalog No. 3247, Catalog 3247I, Catalog No. 3246, and Catalog 3243		
<b>FORMAT:</b>	Liquid concentrate		
<b>PRESENTATION:</b>	One bottle containing 220 mL of PBS, formalin, sucrose and 0.5% sodium azide.		
<b>STORAGE/HANDLING:</b>	Store at 2° to 8°C up to the expiration date printed on the label.		
<b>PREPARATION:</b>	Shake well before each use ensuring that the stock solution is completely mixed. Dilute 1/5 in filtered distilled or deionized water and mix thoroughly.		

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

### *Laboratory Use Reagent*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

© 2002-2006 CHEMICON® International, Inc. - By CHEMICON® International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.