

## HIV - 1 ENV ANTIGEN

<b>CATALOG NUMBER:</b>	AG819
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	100 µg
<b>CONCENTRATION:</b>	1 mg/mL
<b>DESCRIPTION:</b>	HIV-1 ENV (gp160) antigen (recombinant) a.a. 466 to a.a. 753 of the HIV-1 ENV region, 32 kDa.
<b>REACTIVITY:</b>	Reacts strongly with human HIV positive serum.
<b>PURITY:</b>	>95%; Purity of proteins is evaluated by SDS-PAGE, by measuring optical density at 280 nm and by method of Bradford et al.
<b>ORIGIN:</b>	<i>E. coli</i> , recombinant
<b>FUSIONS:</b>	β-galactosidase (114 kDa) fused at the N-terminus
<b>FORMAT:</b>	8 M urea; 20 mM Tris-HCL pH 8.0; 10mM β-mercaptoethanol
<b>STORAGE/HANDLING:</b>	Refrigerate upon arrival. Store at 2 to 8°C for up to 6 months from date of receipt.

*For research use only; not for use as a diagnostic.*

Chemicon offers new Tube-O-DIALYZER sets for dialysis of samples as small as one microliter, with 100% sample recovery. For more information about Tube-O-DIALYZER sets, contact your Chemicon Customer Service Representative today!

Chemicon offers new Tube-O-DIALYZER sets for dialysis of samples as small as one microliter, with 100% sample recovery. For more information about Tube-O-DIALYZER sets, contact your Chemicon Customer Service Representative today!