

**COCAINE- AND AMPHETAMINE RELATED TRANSCRIPT
(CART/CARP)
CONTROL PEPTIDE**

CATALOG NUMBER:	AG228
LOT NUMBER:	
QUANTITY:	50 µg
CONCENTRATION:	1 mg/mL
DESCRIPTION:	Purified mouse/rat CART peptide (55-102).
APPLICATIONS:	For use in blocking the reactivity of AB5338P and AB5340P. Optimal working dilutions must be determined by the end user.
FORMAT:	Purified peptide.
PRESENTATION:	Liquid in PBS, pH 7.5, containing 0.05% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
RELATED REFERENCES:	1. Douglass, J, et al. (1995) <i>J. Neurosci.</i> 15 :2471. 2. Douglass, J and Daoud S (1996) <i>Gene</i> 169 :241. 3. Kristensen, P., et al. (1998) <i>Nature</i> 393 :72.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.