

QUALITY CONTROL RANGES

Milliplex TGF- β 1,2,3 Plex Catalog # TGFB-64K-03

Control Catalog # TGF-6064-3

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	
TGF β 1	Control I	TGFB-3 101	141 - 292	pg/mL
	Control II	TGFB-3 201	546 - 1134	pg/mL
TGF β 2	Control I	TGFB-3 101	146 - 303	pg/mL
	Control II	TGFB-3 201	558 - 1158	pg/mL
TGF β 3	Control I	TGFB-3 101	136 - 282	pg/mL
	Control II	TGFB-3 201	535 - 1110	pg/mL

Millipore • (636) 441-8400 • www.millipore.com