

## ANALYSIS SHEET

### CANINE LEPTIN ELISA CATALOG # EZCL-31K

### CANINE LEPTIN, LYOPHILIZED CATALOG # E8031-K

STANDARD (Lyophilized) LOT No. 831-1K	CANINE LEPTIN VALUE 50 ng/mL
--	---------------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6031-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

### QUALITY CONTROL VALUES CATALOG # E6031-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	CL-101	1.8 – 3.6 ng/mL
QC 2	CL-201	11.3 – 23.4 ng/mL

Millipore  
(866) 441-8400 [www.millipore.com](http://www.millipore.com)

E6031-K.CL.101.201.07.31.08