

## ANALYSIS SHEET

### RAT / MOUSE GHRELIN (ACTIVE) ELISA KIT CATALOG # EZRGRA-90K

#### RAT/MOUSE GHRELIN (ACTIVE) ELISA STANDARD, LYOPHILIZED CATALOG # E8090-K

STANDARD (Lyophilized)	ACTIVE GHRELIN VALUE
LOT No. RGA90-1K	2000 pg/mL*

\* Hydrate standard with 1.8 mL deionized water to obtain 2000 pg/mL value.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6090-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

#### QUALITY CONTROL VALUES CATALOG #E6090-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	ARGR 101	31 – 65 pg/ml
QC 2	ARGR 201	256 - 533 pg/ml

Millipore • (636) 441-8400 • (636) 441-8050 (FAX)  
[www.millipore.com](http://www.millipore.com)

E6090-K.101.201.01.28.08