

ANALYSIS SHEET

BOVINE MORPHICEPTIN RIA KIT CATALOG # 8043K

BOVINE MORPHICEPTIN STANDARD, LYOPHILIZED CATALOG # 8043K

STANDARD (Lyophilized)	BOVINE MORPHICEPTIN VALUE
LOT No. 01BM43	1000 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **1 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # 6043K

Reconstitute each vial with **1 mL** distilled or deionized water.

QUALITY CONTROL VALUES CATALOG #6043K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	BM 101	23.4 – 48.7 pg/mL
QC 2	BM 201	109.0 – 226.3 pg/mL

Millipore
(866) 441-8400 www.millipore.com