

ANALYSIS SHEET

RAT PANCREATIC POLYPEPTIDE RIA KIT CATALOG # RPP-41HK

RAT PANCREATIC POLYPEPTIDE STANDARD, LYOPHILIZED CATALOG # 8041-K

STANDARD (Lyophilized) LOT No. 01P41	RAT PANCREATIC POLYPEPTIDE VALUE 10 ng/mL
---	--

The exact value of lyophilized standard is shown above.

Reconstitute with **1 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # 6041-K

Reconstitute each vial with **1 mL** distilled or deionized water.

QUALITY CONTROL VALUES CATALOG #6041-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	RPP 101	0.5 – 1.1 ng/mL
QC 2	RPP 201	1.5 – 3.0 ng/mL

Millipore
(866) 441-8400 www.millipore.com