

QUALITY CONTROL RANGES

MOUSE SERUM AMYLOID A-3 ELISA CATALOG # EZMSAA3-12K

STANDARD, LYOPHILIZED CATALOG # E8012-K

STANDARD LOT NUMBER: 812-1K	MOUSE SAA-3 VALUE 5.0 µg/mL
--------------------------------	--------------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **0.25 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6012-K

Reconstitute each vial with **0.25 mL** distilled or deionized water.

CONTROL CATALOG # E6012-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	MSAA3 101	0.11 - 0.22 µg/mL
QC 2	MSAA3 201	0.88 – 1.83 µg/mL

Millipore (866) 441-8400 www.millipore.com