

QUALITY CONTROL RANGES

MOUSE ADIPSIN ELISA KIT CATALOG # EZMADPSN-22K

STANDARD, LYOPHILIZED CATALOG#E8022-K

| | |
|--|----------------------------------|
| STANDARD (Lyophilized) LOT No. 822-1K | MOUSE ADIPSIN VALUE 2.5 µg/mL |
|--|----------------------------------|

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6022-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

| CONTROL | LOT NUMBER | EXPECTED RANGE |
|---------|------------|-------------------|
| QC 1 | MADPS 101 | 0.06 - 0.12 µg/mL |
| QC 2 | MADPS 201 | 0.47 – 0.97 µg/mL |

Millipore • (636) 441-8400 • www.millipore.com