

QUALITY CONTROL VALUES
LINCoplex Kit Catalogue # MSCR-42K
Control Catalogue # MSCR6042

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	
sCD30	Control 1	101	605 - 1256	pg/ml
	Control 2	201	6219 - 12916	pg/ml
sgp130	Control 1	101	552 - 1146	pg/ml
	Control 2	201	5445 - 11308	pg/ml
sIL-1RI	Control 1	101	309 - 641	pg/ml
	Control 2	201	2768 - 5749	pg/ml
sIL-1RII	Control 1	101	286 - 595	pg/ml
	Control 2	201	3016 - 6264	pg/ml
sIL-2R α	Control 1	101	269 - 559	pg/ml
	Control 2	201	2643 - 5490	pg/ml
sIL-4R	Control 1	101	294 - 610	pg/ml
	Control 2	201	3135 - 6510	pg/ml
sIL-6R	Control 1	101	560 - 1164	pg/ml
	Control 2	201	6147 - 12768	pg/ml
sRAGE	Control 1	101	549 - 1141	pg/ml
	Control 2	201	5603 - 11637	pg/ml
sTNFR1	Control 1	101	114 - 237	pg/ml
	Control 2	201	1188 - 2468	pg/ml
sTNFR2	Control 1	101	293 - 608	pg/ml
	Control 2	201	2896 - 6015	pg/ml
sVEGFR1	Control 1	101	528 - 1097	pg/ml
	Control 2	201	5911 - 12276	pg/ml
sVEGFR2	Control 1	101	579 - 1202	pg/ml
	Control 2	201	5686 - 11810	pg/ml
sVEGFR3	Control 1	101	577 - 1199	pg/ml
	Control 2	201	5963 - 12384	pg/ml

Millipore • (636) 441-8400 • (636) 441-8050 (FAX)

www.millipore.com

QUALITY CONTROL VALUES
LINCOplex Kit Catalogue # MSCR-42K-1EGFR
Control Catalogue # MSCR6042

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	
sEGFR	Control I	101	559 - 1161	pg/ml
sEGFR	Control II	201	6422 - 13338	pg/ml

Millipore Corporation • (636) 441-8400 • (636) 441-8050 (FAX)
www.millipore.com

MSCR6042.101.201.12.21.07