

ANALYSIS SHEET

RAT/MOUSE OBESTATIN RIA KIT CATALOG # RMOBS-55HK

RAT/MOUSE OBESTATIN STANDARD, LYOPHILIZED CATALOG # 8055-K

STANDARD (Lyophilized)	Rat/Mouse Obestatin VALUE
LOT No. 01RM55	10 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **1 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # 6055-K

Reconstitute each vial with **1 mL** distilled or deionized water.

QUALITY CONTROL VALUES CATALOG #6055-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	OBS101	0.17 – 0.35 ng/mL
QC 2	OBS201	1.04 – 2.16 ng/mL

Millipore
(866) 441-8400 www.millipore.com

6055-K.101.201.04.23.08