

## ANALYSIS SHEET

### CANINE ADIPONECTIN ELISA CATALOG # EZCADP63K

#### CANINE ADIPONECTIN, LYOPHILIZED CATALOG # E8063-K

STANDARD (Lyophilized) LOT No. CAE-1K	CANINE ADIPONECTIN VALUE 100 ng/mL
--	---------------------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6363-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

#### QUALITY CONTROL VALUES CATALOG # E6363-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	CAE-101	0.9 - 1.8 ng/mL
QC 2	CAE-201	7.3 – 15.2 ng/mL

Millipore  
(866) 441-8400 [www.millipore.com](http://www.millipore.com)

E6363-K.101.201.06.07.07