

ANALYSIS SHEET

HUMAN RETINOL BINDING PROTEIN 4 ELISA CATALOG # EZHRBP4-18K

HUMAN RETINOL BINDING PROTEIN 4, LYOPHILIZED CATALOG # E8018-K

STANDARD (Lyophilized) LOT No. 1454972	CANINE ADIPONECTIN VALUE 100 ng/mL
---	---------------------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6018-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

QUALITY CONTROL VALUES CATALOG # E6018-K

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	1454972QC1	4.4 – 9.1 ng/mL
QC 2	1454972QC2	12.6 – 26.2 ng/mL

Millipore
(866) 441-8400 www.millipore.com

E6018-K.1454972.05.21.08