

ANALYSIS SHEET

HUMAN APO AV ELISA CATALOG # EZHAP0AV71K

HUMAN APO AV STANDARD, LYOPHILIZED

STANDARD (Lyophilized)	Human APO AV VALUE
LOT No. E71-1K	370 ng/ml

The exact value of lyophilized standard is shown above.

Reconstitute with **1 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6071

Reconstitute each vial with **1 mL** distilled or deionized water.

QUALITY CONTROL VALUES CATALOG #E6071

CONTROL	LOT NUMBER	EXPECTED RANGE
QC 1	AAV 101	10 – 21 ng/ml
QC 2	AAV 201	54 – 112 ng/ml

Millipore • (636) 441-8400 • (636) 441-8050 (FAX)
www.millipore.com

E6071-K.101.201.04.27.07