

Prostak™ UF Module

Catalogue Number: SK2P001W7

Lot Number: C9DN94311

Membrane Type: Regenerated Cellulose

Membrane Area: 10 ft²

Good Manufacturing Practices

This product was manufactured in a Millipore facility which adheres to Good Manufacturing Practices.

ISO® 9001 Quality Standard

This product was manufactured in a Millipore facility whose Quality Management System is approved by an accredited registering body to the appropriate ISO 9001 Quality Systems Standard.

100% Integrity Testing in Manufacturing

Each Prostak UF module must pass the in-process membrane integrity and housing integrity test.

Quality Assurance Lot Release Criteria

Each module was tested by Manufacturing and released by Quality Assurance to the following specifications:

Integrity

Each module exhibited air flow through fully-wetted membranes less than or equal to 36 cc/min at 10 psi (0.7 bar) inlet pressure at 23 °C.

Housing Integrity

The module housing integrity was tested at 80 psi (5.4 bar) and found to be integral.



Peter Eicher
Quality Manager

Millipore is a registered trademark of Millipore Corporation.
Prostak is a trademark of Millipore Corporation.
ISO is a registered trademark of The International Organization for Standardization.

00107878PU

Rev. B 03/09

